

<< THERMAL METHODS OF ANALYSIS >>


It can be defined as the technique in which change in physical or chemical property of a substance is measured as a function of temperature.

- i> Thermogravimetry or Thermogravimetric Analysis (TG/TGA)
- ii> Differential Thermal Analysis (DTA)
- iii> Differential Scanning Calorimetry (DSC).

* THERMOGRAVIMETRY (TG) OR THERMOGRAVIMETRIC ANALYSIS (TGA)

It is a technique in which change in weight of a substance is measured as a function of temperature or time.

The basic requirements of a thermogravimetry is a precision balance and furnace Assembly for linear rise of temp. with time. Modern thermobalances are computer-controlled and can attain almost any temperature profile. Samples are placed in shallow platinum crucible (sample container) which is connected to a automatic reading microbalance. The most usual type of balance-system is null-point balance system. As there is any change in weight of sample, the needle deflects from its usual position. This deflection is detected by a sensor and a force acts which pulls the needle back to its usual position. The restoring force is proportional to change in weight.







main components of thermobalance

The TG curve for $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ is shown in fig.

i> Horizontal portion signifies regions of no weight loss.

ii> Curve portion signifies regions of weight loss.

The four decomposition rxn are:



The weight loss is due to decomposition of sample at higher temperature.
↳ loss.

The DTG curve is shown in fig 4.

* The region of no-weight loss signifies $dw/dT = 0$.

* The peak in DTG curve represents maximum slope points on the curve.

* In regions B and C, there is change in slope of weight loss curve, these are called inflexion points.


* Inflexion at B arises due to formation of $\text{CuSO}_4 \cdot 3\text{H}_2\text{O}$

* Inflexion at C arises due to formation of $2\text{CuO} \cdot \text{SO}_3$.

* FACTORS INFLUENCING THERMOGRAVIMETRY RESULTS.

i> Heating rate: When a substance is heated at a faster rate, the temperature of decomposition will be higher than that obtained at a slower rate of heating. For a single-step rxn, decomposition curve is shown.

AB = slower rate curve
CD = faster rate curve



$$\begin{aligned} T_A &< T_C \\ T_B &< T_D \\ (T_A - T_B) &< (T_D - T_C) \end{aligned}$$

where T_A & T_C are Temperatures at starting of decomposition.

* Heating rate has less effect for fast reversible rxn.

* The detection of intermediate compounds by TG is very dependent upon heating rate employed.

ii> Furnace Temperature : The nature of atmosphere around the sample [Furnace Atmosphere] have a huge impact on decomposition temperature.

eg: The decomposition of CaCO_3 occurs at much higher rate if carbon is used as atmosphere rather than nitrogen.

The function of atmosphere is to remove the gas products obtained so that the atmosphere throughout the experiment remains as constant as possible. Mostly, vacuum is used.

* Three types of mostly used atmospheres:

i> Static air : Air from surrounding flows through furnace.

ii> Dynamic air : Compressed air from cylinders is passed through furnace.

iii> Nitrogen gas : Oxygen free-nitrogen gas provides inert atmosphere.

NOTE : Sample container should be made up of inert material. eg Platinum, Al.

ii> Crucible Geometry : The shape of crucible can vary the slope of TGA curve. 'flat-plate shaped' are preferred as diffusion of gases evolved is easier.

iv> Sample Characteristics : The weight, particle size of mode of preparation, all alters the TGA results. A large volume of sample can impede (hinder / obstruct) the diffusion of gases evolved. It is better to have small volume of sample and a small particle size.


NOTE : Naturally occurring $\text{Mg}(\text{OH})_2$ has diff. decomposition Temperature than precipitated $\text{Mg}(\text{OH})_2$.

* APPLICATIONS OF ED TGA :

i> determination of purity and chemical stability of primary and secondary standards.

ii> Investigating drying temperatures: Many primary standards absorb water, when exposed in moist air and hence using TGA technique, we can find extent of absorption from which most suitable drying Temperature can be determined.


iii> in qualitative analysis: Any thermogram is a characteristic for a given. Fig. shows TG curve of $\text{CaC}_2\text{O}_4 \cdot \text{H}_2\text{O}$ heated at a rate of $5^\circ\text{C}/\text{min}$.



It is thermally stable upto 1000°C and loses water beyond 100°C . So, CaC_2O_4 is formed which is stable in range $200^\circ\text{C} \rightarrow 346^\circ\text{C}$. At slightly above 346°C , it starts losing CO and CaCO_3 is formed which is stable in range $420^\circ\text{C} \rightarrow 660^\circ\text{C}$. Beyond 660°C , CO_2 is released and CaO is formed which is stable within range $840^\circ\text{C} \& 920^\circ\text{C}$ only.


iv> in study of polymers: It has great application as it gives valuable information about decomposition mechanism of polymers. As polymers have unique pattern, so, it is also used for identification purpose. It is used in examining thermal stability of polymers. Using DTG curve, max. weight change temp. can be determined which is used to distinguish between them.

v> in study of oxidation of alloys:



oxidation of Co_5Sm by TGA in air atmosphere is carried & following curve was obtained. AB represents oxidation of $\text{Sm} \rightarrow \text{Sm}_2\text{O}_3$ & BC represents oxⁿ of $\text{Co} \rightarrow \text{Co}_3\text{O}_4$ and mix compound formation CoSmO_3 . DE represents conversion of $\text{Co}_3\text{O}_4 \rightarrow \text{CoO}$.

The sample container is placed in quartz or pyrex glass within furnace. A thermocouple (A thermoelectric device used to measure temperature & it consists of two metal wires) is located below sample container. The resulting signal is connected directly to x-axis of recorder.



• Block diagram of TGA.

* Characteristics of Good Thermobalance designs.

- i> It should be capable of continuously registering weight change of sample.
- ii> It should reach maximum desirable Temp. ($-150 \rightarrow 2400^\circ\text{C}$ can be obtained with modern thermobalances).
- iii> The rate of heating is linear and reproducible.
- iv> The sample holder should be in hot zone of furnace and this hot zone should have uniform temperature.
- v> The temperature of sample must be measured accurately.

* Methods of expressing TG Results:

The results can be expressed in form of TG curve or DTG curve.

