


Que 1 Explain back titration by taking the example of Volhard's method for determination of chloride ions.

Ans. In chemistry, back titration is a technique used to determine the strength of an analytic through the addition of a known molar concentration of excess reagent. A titration is then performed on the remaining amt. of known solution to determine how much is in excess and to measure the quantity consumed by the analytic.


Volhard's method for determination of chloride ions uses a back titration with KSCN to determine the concentration of chloride ions in a solution. Before the titration an excess chlorine of AgNO_3 soln is added to solution containing Cl^- ions, forming a precipitate of AgCl .



The indicator Fe^{3+} (ferric ion) is then added and the solution is titrated with KSCN solution. The titration remains pale yellow as the excess silver ions react with thiocyanate ions to form silver thiocyanate ppt.



Once all the silver ions have reacted, the slightest excess of SCN^- ions react with Fe^{3+} to form a dark red complex.



The conc. of chloride ions is determined by subtracting the titration findings of the moles of silver ions that reacted with the thiocyanate from the total moles of AgNO_3 added.

Q.2 The percentage transmittance of an solution of unknown compound is 20%, at 25°C & 300 nm for $4 \times 10^{-5}\text{M}$ soln. in a 2 cm long cell. Using this data, calculate

- The absorbance of soln.
- The molar extinction coefficient of compound.

Ans Absorbance of soln = $2 - \log(\%T)$
 $= 2 - \log(20)$
 $= 0.698.$


Using Beer's Law.

$$0.698 = eLC$$

$$= e \times 2 \times 4 \times 10^{-5}$$

$$\Rightarrow e = \frac{0.698 \times 10^5}{8} = 0.087 \times 10^5 \text{ L/mol/cm}.$$

Q.3 The formaldehyde content of a pesticide preparation was determined by weighing 0.3124g of liquid sample into a flask containing 50 ml of 0.0996 M NaOH & 50 ml of 3% H_2O_2 . On heating, the following took place:-



After cooling the excess base was titrated with 23.3 ml of 0.05250 M H_2SO_4 . Calculate the % of HCHO in sample.

Ans Moles of NaOH = $\frac{0.0996 \times 50}{1000}$ or 4.98 mmoles.

Weight of H_2O_2 added = $\frac{3 \times 50}{100} = 1.5\text{g}.$

Moles of H_2O_2 = $\frac{1.5}{34} = 44.1 \text{ mmoles}.$

$$\begin{aligned}\text{Mols of } H_2SO_4 &= 23.3 \times 0.05250 \text{ mmols} \\ &= 12.23 \text{ mmols}.\end{aligned}$$

$$\begin{aligned}\text{Equivalents of } H_2SO_4 (H^+) &= 12.23 \text{ meq} \times 2 \\ &= 24.46 \text{ meq}.\end{aligned}$$

$$\begin{aligned}\text{reacted base} &= (4.98 - 24.46) \text{ meq} \\ &= 2.534 \text{ meq}.\end{aligned}$$

$$\begin{aligned}\text{Hence mols of HCHO} &= 2.534 \text{ mmols} \\ &= 76.02 \text{ mg} = 0.076 \text{ g}.\end{aligned}$$

$$\begin{aligned}\% \text{ strength of HCHO} &= \frac{0.076}{0.3124} \times 100 \\ &= 24.33\%.\end{aligned}$$

Q.4 Describe the principle of differential thermal analysis (DTA) & its comparison with DSC.

Ans Differential Thermal Analysis is a thermanalytic technique. In DTA, the material under study & an inert reference are made to undergo identical thermal cycles, while recording any temp. difference b/w sample & reference. This differential temp. is then plotted against time. Changes in the sample either exothermic or endothermic, can be detected relative to the inert reference. Thus a DTA curve provides data on the transformations that have occurred, such as glass transitions, crystallisation, melting and sublimation. The area under a DTA peak is the enthalpy change and is not affected by the heat capacity of the sample.

DTA

1) The temp. of sample is compared with that of a reference material as both are heated at uniform rate.

2) A DTA output plots $\Delta T = T_R - T_S$ vs temp. of furnace (T_0)

3) It provides calorimetric accuracy of temp. range of -190°C to 1600°C .

4) The area under peak in ordinary DTA is complex fn of sample geometry, heat capacity & heat losses.

DSC


1) The heat energy is supplied at a varying rate to sample or reference, so as to keep their temperatures equal.

2) A DSC output plots heat energy supplied to sample vs identical temp. of two.

3) It provides calorimetric accuracy in temp. range of -170°C to 750°C .

4) The area under a peak can be directly related to enthalpy change occurring.

Q-5 Predict the λ_{max} for following compounds using the Woodward-Fieser rule.




$$\begin{aligned}
 \text{Base value} &= 215 \\
 \text{Extended conjugation} &= 30 \\
 \text{Ring Residue (4)} &= 20 \\
 \text{Exo to Ring (2)} &= 10 \\
 \hline
 &= 275 \text{ nm}
 \end{aligned}$$

20 :-

impl
upli

les.



Ans Base value = 253
 Extended conjugation (2) = 60
 Alkyl substitution (4) = 20
 $\epsilon_{\pi\pi^*}$ (2) = 10
 Ring Residue (4) = $\frac{20}{363 \text{ nm}}$

Q.6 a) Arrange the expected electronic transition for 2-pentanone in order of inc. energy.



2-pentanone is a saturated ketone so $\sigma \rightarrow \sigma^*$ transition must be present ($< 160 \text{ nm}$)

Also $\pi \rightarrow \pi^*$ present, so $n \rightarrow \pi^*$ (189 nm)
 $\pi \rightarrow \pi^*$ (249 nm)


Order of inc energy: $n \rightarrow \pi^* < \pi \rightarrow \pi^* < \sigma \rightarrow \sigma^*$

b) Arrange the following compounds in order of their inc. wavelength of UV absorption

a) Ethylene b) Naphthalene c) Anthracene d) 1,3-Butadiene



Anthracene :



1,3 Butadiene :




Order :- Ethylene < 1,3 Butadiene < Naphthalene < Anthracene.

Q. 7 Indicate which of following vibrations will be IR active or inactive.


Molecule

vibration.



C=C stretching

IR inactive as $\mu = 0$ & $\frac{d\mu}{dn} = 0$.



C-H stretching

IR active as $\mu \neq 0$



N triple N stretching.

IR inactive as $\mu = 0$



Symmetric stretching

The symmetric stretching mode involves no dipole change, so is inactive in IR region.

Q. 8 A solution of sodium hydroxide contained 0.250 mol cm^{-3} . Using phenolphthalein indicator titration of 25 cm^3 of this solution required 22.5 cm^3 of a hydrochloric acid solⁿ for complete neutralisation.

a) Write the eqⁿ for titration.



b) what apparatus would you use to measure out

i) the sodium hydroxide solⁿ.

Ans. volumetric pipette.

ii) HCl solⁿ.

Ans buret.

c) what would you rinse your apparatus but with before doing the titration.

Ans water

d) Calculate the moles of NaOH neutralised


$$\begin{aligned}\text{No. of moles of NaOH} &= 0.250 \times 25.0 \\ &= 6.25 \text{ mmols.}\end{aligned}$$

e) Calculate the conc. of HCl in mol/dm³.


$$\text{conc. of HCl} = \frac{6.25 \times 10^{-3}}{22.5} = 0.277 \text{ mol/dm}^3$$

Q.9 At pH > 12, phenolphthalein becomes colorless. Explain with the help of molecular structure.

Ans In mild basic conditions, the structure of phenolphthalein is



But when excess of OH^- ions are present.




Q.10 In iodometric titration, what should be nature of analyte? Explain with chemical rxns.

Ans In iodometric titration, the nature of analyte must be oxidative.


To a known volume of sample, an excess but known amt. of iodide is added, which the oxidizing agent oxidises iodide to iodine.



The iodine dissolves in the iodide-containing solⁿ to give triiodide ions, which have a dark brown colour.



The triiodide ion solⁿ is then titrated against std. thiosulphate to give iodide again.



The overall rxn is thus.

